

KT-TD3 Portable Bench Top Low Speed Centrifuge (PRF CGF dental centrifuge)



General Introduction

KT-TD3 beauty centrifuge is specialized in the purification of FAT and PRF / PRP / CGF, in order to improve the survival rate of FAT, a lot of research has been done in terms of centrifugal speed, time, centrifugal force, diameter, etc. The multi-functional purification centrifuge improves the efficiency of surgery for autologous blood skin regeneration, shortens the operation time, maximizes the survival rate of lipid and PRP during the operation, makes the transplant operation simple and convenient, and is a good choice for plastic surgeons.

Application

beauty clinics, cosmetic surgery hospitals, beauty salons, health centers



Technical Data

Model	KT-TD3
Max Speed	4000 rpm
MAX RCF	2200xg
MAX Capacity	8*15ml
Speed Accuracy	±20rpm
Timer Range	0~99min
Display	LED
Motor	DC brushless
Power	50W
Power Supply	AC 220V 50HZ 2A
Net Weight	6KG
Device Dimension	320*260*230(mm)
Packing Size	340*280*250(mm)

Matched Rotor

NO	ITEM	MAX SPEED(rpm)	MAX RCF(xg)	Capacity(ml)	Tube Size(mm)
NO.01	Angle rotor	4000rpm	2280xg	6*10/15ml	Ø 16*120
NO.02	Angle rotor	4000rpm	2280xg	8*10/15ml	Ø 16*120
NO.03	Angle rotor	4000rpm	2280xg	12*5/7ml	Ø 13*95
NO.04	Angle rotor	4000rpm	2280xg	6*20ml	Ø 20*120