


Inoculating Loop & Needle Sterilizer KT-IS800B (ϕ 35mm)



Description:

KT-IS800B Inoculating Loop & Needle Sterilizer is used in degassing for inoculating loops and needles, instead of alcohol lamps. Due to its convenience, sterilizer is widely used in the environment of biosafety cabinets, clean benches. They are also used besides exhaust fans or on mobile vans.

Features:

- ◆ Safety, convenient. Totally take the place of alcohol lamps.
- ◆ Only 5 seconds to 7 seconds are needed, because of its highest temperature of 825°C.
- ◆ Smart, light, easy to handle.
- ◆ Used in anaerobic chamber is allowed.
- ◆ Podzolic organics prevent pollutions effectively.
- ◆ Smart technique makes it possible to use for a long time and much safe.

Application:

Widely used in biosafety cabinets, clean bench, near exhaust fan, mobile vehicles, even outdoor and other severe environment can be easily used.

Parameters:

Parameter	Type	KT-IS800B
Max. temp. in central zone		825°C±25°C
Max. outer diameter of sterilization goods		ϕ 35mm
Length of heating zone		100mm
Heating time		15min
Ambient temperature		5°C~40°C
Fuse		250V 3A ϕ 5×20
Voltage		AC220V or AC110V, 50-60Hz
Power		250W
Dimension(mm)		130×120×200
Weight(kg)		1.3