

KT-TXK4 Blood Type Card /Gel Card Centrifuge



Features

1. Small size, save laboratory space.
2. Adopt all steel structure, stainless steel centrifugal cavity.
3. Maintenance-free DC brushless motor drive, smooth and quiet operation.
4. The parameters can be changed at any time during operation without stopping.
5. Automatic calculation and synchronous display of centrifugal force RCF value.
6. The use of electronic door locks enhances safety.
7. This model is widely used in laboratories at all levels of hospitals, colleges and universities due to its high cost performance and stable and reliable quality.

Technical Data

Model	KT-TXK4
Max Speed	4000 rpm
Max RCF	2240 xg
Max Capacity	12/24 cards
Speed Accuracy	± 50rpm
Timer Range	1 - 99min
Power	200W AC 220V 50HZ 10A
Noise	≤55dB(A)
Net Weight	24Kg
Device Dimension	483*320*265mm (L*W*H)